# Wrist Rehabilitation

## cold compression therapy

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**Contents**  
Nylon, Polyester, Polycarbonate, Polypropylene, Polyethylene, Polyvinyl chloride  
**Gelpack**  
Propylene glycol, Sodium carboxymethyl cellulose

Read the instructions for use before using the product. Lees voor gebruik van dit medisch hulpmiddel de gebruiksaanwijzing. Veuillez lire attentivement la notice d'utilisation avant d'utiliser ce produit. Lesen Sie die Gebrauchsanweisung sorgfältig durch bevor Sie dieses Produkt verwenden.

**GB** Incorrect application can cause injuries. Consult a qualified person for the correct application. Stop using the product immediately when skin irritation occurs. Keep out of reach of children. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Store in a well ventilated place. Dispose of contents / container in accordance with local / regional / national / international regulations.

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
GB

#### PARTS

- 1 compressor-cover
- 2 gelpack
- 3 pump with tube
- 4 vent button

#### INSTRUCTION FOR USE

1. Put the gelpack (2) in the fridge/freezer at least 2 hours before use.
2. Put the gelpack (2) with its back (velcro tape) on the compressor-cover (1).
3. Place the thumb (of the injured wrist) through the small loop (A) and place it on your leg.
4. Winding around the wrist and pull the flap (image 1) on and over the cover and attach the velcro tape in place (image 2).
5. Connect the tube (3) to the ON/OFF switch, screw it on the two-way valve.
6. Turn the switch in the ON position.
7. Pump it up (for compression) and increase the pressure around the affected wrist\*.
8. Turn the switch to the OFF position and just twist and unwind the tube.
9. Brace on for 5-10 minutes. Brace off for 5-10 minutes and repeat.
10. Empty the compressor-cover before taking the brace off: Push in the vent button (4).

\*Avoid over-compression; excessive pressure may cause blood circulation disorder.

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